

C6 Carriage for SC-DX Autosamplers

Highly Accurate and Durable Carriage for SC-DX Autosamplers

The new and improved SC-DX autosamplers are now available with the C6 carriage, a durable, long-life carriage that is designed with a ridged column and Reset Probe for the highest level of accuracy. The C6 carriage is topped with a newly-designed 360 cap that makes for easy probe installation and a shorter length of capillary.

BENEFITS:

- **Durable and Accurate** - Ridged vertical column facilitates precise movement
- **Reset Probe** - The Reset Probe protects against probe damage and ensures no samples are missed
- **360 Cap** - Easy probe installation and shorter capillary length
- **Faster Sampling** - Light-weight design and smooth vertical motion allows for faster sampling
- **Long Lifetime** - Chemically-resistant construction with Teflon® coating ensures accurate sampling
- **Upgrade Available** - Existing SC-DX autosamplers can be upgraded to the new C6 carriage
- **Septum-Piercing Option Available**
- **Optional Configuration for Vessels up to 190 mm Tall**




C6 Carriage with 360 Cap

The C6 carriage is now standard on the following autosamplers:

- SC-2 DX
- SC-4 DX
- SC-8 DX
- SC-14 DX
- SC-E2 DX
- SC-E4 DX

SC-DX autosamplers with the 6000 carriage can be upgraded to the C6 carriage.



SC-2 DX FAST Autosampler with C6 Carriage

C6 Carriage Upgrade Part Numbers		
Standard (140 mm)	Tall (165 mm)	Description
YZ-C602	YZ-C602-18	SC-2 DX and SC-4 DX Carriage
YZ-C608	YZ-C608-18	SC-8 DX and SC-14 DX Carriage
YZ-C652	YZ-C652-18	SC-E2 DX and SC-E4 DX Carriage for Element 2

Standard: For sample vessels up to 140 mm tall.
 Tall: For sample vessels up to 165 mm
 Z-Axis Extension arm available for 190 mm tall sample vessels.