

## FAST DXi

Integrated FAST Valve and Precision Micro Peripump for ICP/ICPMS

- Direct replacement for standard pump
- Connects to any SC-DX autosampler

### Precision

- Accurate pumping across a wide range of rates from low nL/min to 10 mL/min
- Faster analysis with high precision.

### Low Pulsation

- 12 closely-spaced rollers maintain consistent flow rate

### Chemical Resistance

- Ceramic pins and Teflon®-coated metal components are resistant to acids and organic solvents

### Hardware Advantages

- Four channel pump
- Pump Design Extends Tubing Lifetime
- 1, 2, 3, 4 and 6 Channels

### Peripump Control Options

- Instrument software control or
- SC-FAST software control



### APPLICATIONS:

- Precise On-line Internal Standard Addition



- Low-Flow Applications



- Microchromatography



- High Precision - SSI

