

PFA PureTip Pipette Tips

An essential component for ultratrace analysis.


Often overlooked as a source of sample contamination the pipette tip can come into direct contact with the sample or acids that are added to the sample.

Only PureTips can deliver contamination free concentrated chemicals to the sample. PureTips can be reused hundreds of times and will not leach contaminants throughout their lifetime.

- Molded from high purity PFA fluoropolymer, minimizing contamination
- Chemically resistant to virtually all sample types, ideal for strong acids, alkalis, organics and biologicals
- Long-life, reusable
- Use with standard Eppendorf type 1000 μ L & 100 μ L pipetters

Ordering Information

Part Number	Volume	Qty	Description
ES-7000-0101	100 μ L	10 pk	Pre-cleaned PFA Pure Tip pipette tips 100 μ L
ES-7000-1001	1000 μ L	10 pk	Pre-cleaned PFA Pure Tip pipette tips 1000 μ L



PureTip Pipette Tips are available in 1000 μ L & 100 μ L sizes and can be used with standard Eppendorf type pipetters.